

HROMATEC



DRP 4 Manual Syringe Headspace for GC

DRP4 - economical solution enables you to use the headspace technique to measure volatile and semi-volatile compounds in your sample without additional sample preparation. The system does not connect directly with the GC and it can be used with any existing gas chromatography system.

The high cost of automated headspace systems limited the headspace technique uses, but now, DRP4 with low ownership cost and advanced design archives excellent performance with high accuracy and reproducibility.

Application range:

Manual headspace DRP4 is an ideal tool to analyze compounds with low boiling point temperatures for low and medium throughput laboratories and can be used for environmental monitoring, forensic and pharmaceutical analysis, oil and petrochemicals applications.

- Hydrocarbon traces in soil
- Solvent traces in drugs
- Alcohol in blood
- Flavor analysis in food
- Light Hydrocarbons in oil

- VOC analysis in water
- Methanol/Ethanol in Drinks
- Volatile compounds in plastic & polymers
- Residual Solvent Analysis
- and many others

DRP4 Key features:

- Standalone and compatibly with any gas chromatograph
- Low ownership and service cost
- Syringe purging by inter gas
- Easy to use
- CPU controller with Integrated timer
- Air Removing from the vial

Temperatures Sensors:

The temperature modes in headspace analysis are very important to get reliable results. And DRP4 has a possibility for digital calibration of all temperatures sensors with accuracy: 0.1 °C





She echnic

Technical Specification:

General Specification:

Operating principle Syringe type: Vial incubator working principle: User interface:

Syringe heating:

Volume Sample volume range: Syringe temperature: Syringe temperature setpoint resolution: Syringe cleaning : Syringe purging time: Needle:

Vial Incubator:

Vial incubator capacity: Vial size: Vial incubator temperature: Incubator temperature setpoint resolution: Shaking method: Shaking speed: Shaking cycles : Vial incubation time: Analytical cycle time:

Environmental Conditions:

Ambient Operating Temperature: Relative humidity: Storage Temperature: Power Requirements: Power consumption:

Other specification:

Dimensions: (WxDxH); Weight: Syringe headspace sampling Special headspace SGE Diamond [™] Syringe Heating & Shaking 4-line LCD display, complete monitoring and control parameters

1 or 2.5 ml 0.1 – 2.5 ml 35 ... 150°C 0,1 °C Purging by inert gas from the top through syringe barrel 1 ... 99 min Fixed in the syringe barrel, side hole needle tip

4 vials 20 mL and 10 mL 35 ... 150°C 0,1 °C Orbital Fast / Slow / Off On / Off, 1 ... 999 sec 1 ... 999 min 1 ... 999 min

from 10 to 35°C not more than 80 % from -50 °C to 50°C ~220V ±10%, 50Hz 150 W

245 mm x255 mm x 150 mm 7 kg

Safety and Certification:

Products designed and manufactured under regulations of GOST R ISO 9001 quality standard. At electromagnetic compatibility the chromatograph meets the requirements of IEC 61010-1

Information and technical specification in this publication are subject to change without notice.

© Chromatec JSC SDO, 2015 November 24, 2015 09-101-5012EN



For more information please contact us at: info@chromatec.ru